

# SOMNI 3 VAPORIZER



## LIFE CRITICAL

- ✓ Low flow calibration at 21 degrees centigrade using **LASER REFRACTOMETRY** ensures un-paralleled accuracy of anesthetic vapor delivery throughout the dial percentage range.
- ✓ Brass and specific alloy materials assure accurate **HEAT TRANSFERENCE** which is **CRITICAL** to control anesthetic agent **VAPORIZATION**.
- ✓ The internal **ADVANCED WICKING** system (*currently utilized in human anesthetic vaporizers*) assures linear vapor delivery during **EXTENDED** surgical procedures.

## QUALITY

- ✓ External and internal parts are machined from **BRASS and SPECIFIC ALLOYS**, which are proven **NOT** to react with anesthetic agents, which has been identified to be a corrosion problem with certain aluminum and plastic materials.
- ✓ **STAINLESS STEEL** cagemounts ensure durability over many years and procedures.

## SAFETY

- ✓ The liquid holding sump machined from one solid piece of **BRASS** assures many years of **SAFE LEAK PROOF** performance.
- ✓ The **ONE PIECE FILLING SYSTEM** simplifies the filling process and eases agent level monitoring during procedures.

## SPECIFICATIONS

Weight	15.5 lb / 7 kg
Height	7.58 in / 195 mm
Depth	5.5 in / 140 mm
Width	5.25 in / 135 mm
Agent Capacity	250 ml
Wicking Capacity	75 ml

**SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.**

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625    ✉ info@somniscientific.com

🌐 www.somniscientific.com



*Proudly machined, assembled  
and calibrated in the USA*