



O2 CONCENTRATOR

QUALITY

- ✓ Triple air filtration including carbon and bacterial filtration.
- ✓ Multi parameter display, including percentage of concentration.
- ✓ Easy filter access.
- ✓ Consistent supply of high concentration supplemental oxygen.

ECONOMICAL

- ✓ Has the ability to pay for itself in less than a year.
- ✓ Low total cost of ownership.
- ✓ Eliminates the requirement of delivered oxygen.

SAFE

- ✓ Power failure, no flow, high and low pressure alarms.

CONVENIENT

- ✓ Smooth quiet operation.
- ✓ Casters for portability.
- ✓ Programmable digital run timer.
- ✓ Easily adaptable to any anaesthesia system (O2 NIST MALE and Hose Barb fitting included)

SPECIFICATIONS

Weight	5 lpm 27 kg OR 10 lpm 29 kg
Height	63.5 cm
Depth	31.75 cm
Width	38.1 cm
Flow Meter Capacity	0-5 lpm OR 0-10 lpm
Sound	40 - 45 dBa
Line Pressure	10 psi

(T) 0800 0129101 (D) 01872 248890 (M) 07798 969805
 enquiries@somniscientific.com www.somniscientific.co.uk

INHALATION ANAESTHESIA EQUIPMENT AND SERVICE

SOMNI Scientific is centered on the animal welfare and research community with a focus on clinical accuracy, clinician/technician safety, economic performance and intuitive functionality.

