

Quip Laboratories, Inc. MB-10 Tablets are the most convenient way to deliver the benefits of Chlorine Dioxide.



QUALITY

- ✓ Virucidal, Bactericidal, Tuberculocidal
- ✓ Featured on the Center for Biocide Chemistries' list of Novel Coronavirus (COVID-19)—Fighting Products.

SAFETY

- ✓ Safe and Stable, zero toxicity
- ✓ Relatively neutral pH balance and excellent material compatibility

CONVENIENT

- ✓ Tablet shelf -stable for 3 years and 7 days once dissolved
- ✓ Dissolves within minutes, disinfects within minutes
- ✓ Recommended for Disinfection of SOMNI InductFlow Red and Clear Induction Chambers



SPECIFICATIONS

ONE MB-10 1.5 GRAM TABLET DISSOLVED IN 500 ML OF WATER = 200PPM CHLORINE DIOXIDE DISINFECTANT SOLUTION

ITEM# CC-3110 SUPPLIED AS A PACK OF 10
ITEM# CC-3125 SUPPLIED AS A PACK OF 25

(T) 0800 0129101 (D) 01872 248890 (M) 07798 969805
enquiries@somniscientific.com www.somniscientific.co.uk

INHALATION ANAESTHESIA EQUIPMENT AND SERVICE

SOMNI Scientific is centered on the animal welfare and research community with a focus on clinical accuracy, clinician/technician safety, economic performance and intuitive functionality.

