

Somni Scientifics AMD 3+ (Accurate Multiple Delivery) is the ultimate in Inhalant Anaesthesia gas delivery. This pressurised anaesthesia system can deliver accurate adjustable flowrates to 3 stations simultaneously. Each station is equipped with a flowmeter and toggle switch that allows the user to leave flows preset or adjust to their needs. Unlike any other system on the market this allows users the greatest flexibility for their procedures. Any combination of Nose Cones, Induction Chambers, Stereotaxic or Imaging device can be used with the AMD-3+.



FEATURES

- ✓ 3 Flowmeter Controlled Stations with Toggle Switches
- ✓ Color Coded Tubing to each station
- ✓ Common Gas outlets for each station
- ✓ Small footprint
- ✓ Adjustable Flowrates
- ✓ Stacking design saves valuable workbench space

INCLUDES

- ✓ TT-2 with Special 6psi Calibrated vaporiser
- ✓ 3 each, 6' color coded tubing with adapters

SPECIFICATIONS

Weight	
Height	7 in / mm
Depth	8 in / mm
Width	12 in / mm

SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.

For more information or pricing, please contact our sales team.

☎ (T) 0800 0129101 (D) 01872 248890 (M) 07798 969805

✉ enquiries@somniscientific.com 🌐 www.somniscientific.co.uk