


INDUCTFLOW CHAMBERS

ACTIVE VACUUM AND PASSIVE STYLES

The SOMNI InductFlow line of induction chambers are designed with slide top lids that allow for one handed operation while preventing the anaesthetic gas from being drawn into the workplace. The red tinted acrylic provides for clear subject observation while reducing subject stress during anaesthetic procedure.



QUALITY

- ✓ Red tint creates less stressful environment.
- ✓ Removable floor liner.
- ✓ Constructed of shatter-resistant acrylic.

SAFETY

- ✓ Prevents unnecessary exposure to personnel and work environment (NIOSH 2ppm) when used with active waste gas removal system EPS-3 or W.A.G. filter (active chamber only).

CONVENIENT

- ✓ One handed operation.
- ✓ Color coded inlet and outlet ports.
 - Compatible with 19mm tubing & 15mm adapters.

SPECIFICATIONS

ACTIVE

Medium	9.25 in L	4.25 in W	5.75 in H
Large	11.75 in L	6 in W	5.75 in H

PASSIVE

Small	4.25 in L	4.25 in W	4.25 in H
Medium	9.25 in L	4.25 in W	4.65 in H
Large	11.75 in L	6 in W	5.75 in H

SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.

For more information or pricing, please contact our sales team.

☎ (T) 0800 0129101 (D) 01872 248890 (M) 07798 969805

✉ enquiries@somniscientific.com 🌐 www.somniscientific.co.uk