

SOMNI
SCIENTIFIC



SOMNI BAIN MANIFOLD



SAFETY:

- Pressure Manometer center mount for immediate visibility
- Popoff valve releases excess pressure at 2cm H₂O while maintaining passive volume in the rebreathing bag
- Monitoring pressure during manual ventilation is critical in preventing barotrauma or inadequate ventilation, especially in small patients

Standard anaesthesia systems do not incorporate this feature

► SPECIFICATIONS

WEIGHT: .90 kg.

HEIGHT: 20.32 cm.

DEPTH: 19.05 cm.

WIDTH: 7.62 cm.

SOMNI Bain Nonrebreathing
Assembly shown with RM-17003
Small Animal Rebreathing
Anaesthesia Machine and VS-6001
SOMNI 3 Sevoflurane Vaporiser
Funnel Fill 250 ML Capacity




CONVENIENT:

- Simple press and twist operation to anchor.
- Freely rotates on the base.
- Easy and quick to move when flipping the patient.



QUALITY:

- Manufactured from high quality, durable materials
- Adjustable manometer
- Space conscious
- Easily adapted to most breathing systems



(T) 0800 0129101 (D) 01872 248890 (M) 07798 969805

enquiries@somniscientific.com www.somniscientific.co.uk