

SOMNI
SCIENTIFIC



AMD 3

ACCURATE MULTIPLE DELIVERY

Somni Scientific's AMD (Accurate Multiple Delivery)

The Somni AMD-3 accurate multiple delivery is a flowmeter controlled diverter system that is adaptable to any standard Inhalant Anaesthesia System. It can be added to Somni's TT-2 or any other non-pressurized system.

The AMD-3 gives the user the ability to deliver accurate flowrates to 3 stations simultaneously. No more guessing if the flow is being delivered equally to multiple devices at once. With its included pressure relief valve the user does not have to worry about pressure building up in the vaporiser. 3 fresh gas common outlets supplied with color coded tubing makes this system very user friendly.

► SPECIFICATIONS

Weight: 11 lbs / 5 kg

Height: 7 in / 177.8 mm

Depth: 8 in / 203.2 mm

Width: 12 in / 304.8 mm



FEATURES:

- 3 Flowmeter Controlled Stations with Toggle Switches
- Color Coded Tubing to each station
- Common Gas outlets for each station
- Small footprint
- Adjustable Flowrates
- Stacking design saves valuable workbench space




INCLUDES:

- Pressure Relief Valve
- 3 each, 6' color coded tubing with adapters



(T) 0800 0129101 (D) 01872 248890 (M) 07798 969805

enquiries@somniscientific.com  www.somniscientific.co.uk