

SOMNI
SCIENTIFIC



O2 CONCENTRATOR



QUALITY:

- Triple air filtration including carbon and bacterial filtration.
- Multi parameter display, including percentage of concentration.
- Easy filter access.
- Consistent supply of high concentration supplemental oxygen.



SPECIFICATIONS

WEIGHT: 5 lpm: 35 lbs OR 10 lpm: 38 lbs

HEIGHT: 635 mm / 25 in

DEPTH: 317.5 mm / 12.5 in

WIDTH: 381 mm / 15 in

FLOW METER CAPACITY: 0-5 lpm OR 0-10 lpm

SOUND: 40 - 45 dBa

LINE PRESSURE: 10 psi

POWER OPTIONS: 220V



ECONOMICAL:

- Has the ability to pay for itself in less than a year.
- Low total cost of ownership.
- Eliminates the requirement of delivered oxygen.




SAFE:

- Power failure, no flow, high and low pressure alarms.



CONVENIENT:

- Smooth quiet operation.
- Casters for portability.
- Programmable digital run timer.



(T) 0800 0129101 (D) 01872 248890 (M) 07798 969805

enquiries@somniscientific.com www.somniscientific.co.uk