

SOMNI
SCIENTIFIC



SOMNI C3 125 ML VAPORISER KEY-FILL

SOMNI C3 125 ML VAPORISER KEY-FILL

The SOMNI C3 is a precision medical vaporiser designed for Veterinary Use Only. This vaporiser was designed to mix the vapours of specified liquid anaesthetic gas agents with anaesthetic delivery gases.



QUALITY:

- External and internal parts are machined from BRASS and SPECIFIC ALLOYS, which are proven NOT to react with anaesthetic agents, which has been identified to be a corrosion problem with certain aluminium and plastic materials.
- STAINLESS STEEL cagemounts ensure durability over many years and procedures.

SPECIFICATIONS

WEIGHT: 6.1kg / 13.5 lb

HEIGHT: 181.6 mm / 7.15 in

DEPTH: 134.6 mm / 5.35 in

WIDTH: 109.2 mm / 4.3 in

CAPACITY: 125 ml

WICKING CAPACITY: 50 ml



SAFETY:

- The liquid holding sump machined from one solid piece of BRASS assures many years of SAFE LEAK PROOF performance.
- The ONE PIECE FILLING SYSTEM simplifies the filling process and eases agent level monitoring during procedures.



LIFE CRITICAL:

- Low flow calibration at 21 degrees centigrade using LASER REFRACTOMETRY ensures un-paralleled accuracy of anaesthetic vapor delivery throughout the dial percentage range.
- Brass and specific alloy materials assure accurate HEAT TRANSFERENCE which is CRITICAL to control anaesthetic agent VAPORISATION.
- The internal ADVANCED WICKING system (currently utilised in human anaesthetic vaporisers) assures linear vapour delivery during EXTENDED surgical procedures.



(T) 0800 0129101 (D) 01872 248890 (M) 07798 969805

enquiries@somniscientific.com www.somniscientific.co.uk